

Contrast Imaging Phantom

Model 07-643



- Check dynamic range of the video system
- Check the system all at once
- Use it during all fluoro modes - pulsed, non-pulsed, etc.
- Check film range and density

Introduction

The Contrast Imaging Phantom (Model 07-643) is an accurate, easy-to-use, indispensable tool for evaluating image quality and determining that the imaging system is operating at its full potential. It will immediately let you know if there's a problem.

Specifications

Dimensions


Outside diameter 9.05 in (23 cm)

Thickness 0.5 in (1.28 cm)


Weight 2.80 lb (1.26 kg)

Available model(s)

07-643 Contrast Imaging Phantom



Model 07-643 diagram



Model 07-643 x-ray image

We now have available sterile Image Intensifier Covers, call for details.

For more information, receive our full product catalog, or order online, contact **Radiation Management Services** business of **Fluke Biomedical**: 440.248.9300 or www.flukebiomedical.com/rms.

Specifications are subject to change without notice.

©2005 Fluke Biomedical. All rights reserved. Printed in USA. 07-643-ds rev 2 22 feb 05