



Manipulation tong for remote handling of radioactive sources

Ré f. POSITONG



Positong 350 represented with its Posishield.



Principle

The ergonomic Positong was designed to handle vials and syringe while protecting the operators from undue exposure risks. Its ergonomically designed handle, provides comfortable gripping and integrates an adjustment knob to limit the opening size of the jaws.

As an option, in order to reduce exposure to the hand when handling high energies and/or high activity sources, the Positong can be equipped with a 3 mm thick tungsten 'Posishield'.

This is excellent when handling FDG or I131

Technical description

- Ergonomic plastic handle, stainless steel jaws.
- Easily decontaminable unit.
- Adaptable for right or left handed use.
- Selectable jaw opening screw knob.

• Option

Removable 'Posishield' 3 mm thick tungsten shield to be screwed on the handle.

Ergonomic Positong

Ref. Positong 100 Length (L) 100 mm

Ref. Positong 200 Length (L) 200 mm

Ref. Positong 350 Length (L) 350 mm

Ergonomic shield Posishield

Ref. Posishield

For Positong 200 & Positong 350

Size: L.80 x H.130 mm

