

Acrylic beta syringe shield


Type beta

Principle

The clear acrylic beta syringe shields allow 360° viewing of the syringe contents while providing maximum beta ray protection to the user. The beveled shape and tip gives optimum manoeuvrability during vinipuncture.

Technical description

- Protection by use of acrylic 10 mm thick
- Fast syringe gripping system with push button.



Reference	Syringes	Protection length	Inside Ø	Weight
BETA 1	1 cc	75 mm	9 mm	45 g
BETA 2	2 cc	52 mm	11.5 mm	38 g
BETA 2.5	2.5 cc	58 mm	11.5 mm	45 g
BETA 3	3 cc	67 mm	11.5 mm	51 g
BETA 5	5 cc	65 mm	15 mm	55 g
BETA 10	10 cc	83 mm	18 mm	80 g



3, rue de l'Europe - 44470 CARQUEFOU (France) - Tel. 33 (0)2 40 25 24 04 - Fax 33 (0)2 40 25 18 37
e-mail: contact@lemerpax.com - www.lemerpax.com