

# Extension Cable Reel for Ionization Chambers

## Model 30-355

- Prevents tangled, loose, or misplaced triaxial-connector cables
- A must for every therapy department

The Model 30-355 cable reel offers a great convenience in winding and storage of extension cables. This cable reel will extend the life of a cable significantly by eliminating kinks and providing protection during storage. Only as much cable as needed has to be reeled out at any time.

A cable reel will accommodate a white extension cable up to 50 meters in length or a grey cable up to 20 meters. Either gender may be specified at the hub, depending on whether the cable reel will be kept at the electrometer (male connector at the hub) or inside the treatment room (female connector at the hub).

**Dimensions** 8.5 in Ø x 3.0 in (h) (21.5 x 7.5 cm)

**Weight** 1 to 3 lb (depending on the cable length)

**Available model(s)**

**30-355** Extension Cable Reel for Ionization Chambers



---

For more information or to receive our full product catalog, contact **Fluke Biomedical** at 440.248.9300 or [www.flukebiomedical.com/rms](http://www.flukebiomedical.com/rms).

*Specifications are subject to change without notice.*

©2005, 2006 Fluke Biomedical. All rights reserved. Markus and Roos are trademarks of PTW New York and PTW Freiburg. Printed in USA.  
30-331-ds rev 3 06 mar 06