

# Lead Aprons

Models 57-940,  
57-941 and 57-959 Series

Available are three highly popular lead apron styles to meet your shielding and budget needs: the Thyro-Shield Lead Apron, the Tie-Closure Lead Apron, and the Belt-and-Buckle Lead Apron. All three styles feature shoulder and back supports for added comfort. They are fabricated from high-quality flexible materials that allow them to drape easily, so they don't restrict body movement. All aprons are available in dark blue, light blue, and beige.



- Comfortable, effective protection for personnel working in high radiation areas
- Maximum radiation protection in one convenient apron
- Nonrestrictive design
- Nonabsorbent and easy to clean
- Resists chemicals, abrasion, and aging
- Three styles and colors to choose from

## Thyro-Shield™ Lead Apron

## Tie-Closure Lead Apron

Model	Lead equiv.	Width (inches)	Length* (inches)	Weight
57-959-3218	0.5 mm	18	32	10.5 lb (4.8 kg)
57-959-3400	0.5 mm	24	34	11.7 lb (5.4 kg)
57-959-3600	0.5 mm	24	36	12.5 lb (5.7 kg)
57-959	0.5 mm	24	38	13.5 lb (6.2 kg)
Model numbers shown are for dark blue only. To order light blue or beige, request model number.				
* Length is from neckline to bottom of apron.				

Model	Lead equiv.	Width (inches)	Length* (inches)	Weight
57-941-3218	0.5 mm	18	32	9.7 lb (4.4 kg)
57-941-3200	0.5 mm	24	32	10.2 lb (4.6 kg)
57-941-3400	0.5 mm	24	34	10.4 lb (4.7 kg)
57-941-3600	0.5 mm	24	36	11.1 lb (5 kg)
57-941	0.5 mm	24	38	12.2 lb (5.4 kg)
Model numbers shown are for dark blue only. To order light blue or beige, request model number.				
* Length is from neckline to bottom of apron.				

## Belt-and-Buckle Lead Apron

Model	Lead equiv.	Width (inches)	Length* (inches)	Weight
57-940-3218	0.5 mm	18	32	9.7 lb (4.4 kg)
57-940-3200	0.5 mm	24	32	10.2 lb (4.6 kg)
57-940-3400	0.5 mm	24 <td 34	10.4 lb (4.7 kg)	
57-940-3600	0.5 mm	24	36	11.1 lb (5 kg)
57-940	0.5 mm	24	38	12.2 lb (5.4 kg)
Model numbers shown are for dark blue only. To order light blue or beige, request model number.				
* Length is from neckline to bottom of apron.				



## Thyro-Shield™

Model 57-753

*Reduce radiation exposure to the thyroid*

- Lightweight, comfortable, easy to use
- Shields the thyroid area with 0.5 mm lead-equivalent protection

**Dimensions** 23 in (l), 2 to 4.75 in (w)  
**Weight** 0.43 lb (0.2 kg)



For more information, receive our full product catalog, or order online, contact **Radiation Management Services** business of **Fluke Biomedical**: 440.248.9300 or [www.flukebiomedical.com/rms](http://www.flukebiomedical.com/rms).

*Specifications are subject to change without notice.*

©2005 Fluke Biomedical. All rights reserved. Thyro-Shield is a trademark of Fluke Corporation. Printed in USA.  
57-940-ds rev 2 10 jun 05