

Small shielded safes

Principle

The small shielded safes are designed to the storage of radioactive sources as syringes, capsules, vials and standards. They are equipped with a key lock.

Technical description

- White decontaminable finish.
- Key locked latch.
- Rubber seal on the opening.



Reference	Lead thickness	Inside dimensions	Outside dimensions	Weight
SS25/200X200	25 mm	L.200 x H.200 x D.200 mm	L.313 x H.320 x D.265 mm	120 kg
SS25/300X300	25 mm	L.300 x H.300 x D.300 mm	L.413 x H.420 x D.365 mm	210 kg
SS25/300X500	25 mm	L.300 x H.500 x D.300 mm	L.400 x H.640 x D.370 mm	325 kg
SS40/300X400	40 mm	L.300 x H.300 x D.400 mm	L.445 x H.530 x D.410 mm	470 kg
SS40/500X575	40 mm	L.500 x H.600 x D.575 mm	L.600 x H.750 x D.600 mm	1100 kg

Other sizes on request.